

LENS PARAMETERS

MATERIAL:	samfilcon A
LENS MATERIAL TECHNOLOGY:	MoistureSeal®
WATER CONTENT:	46%
OXYGEN TRANSMISSION:	163 Dk/t @ center for -3.00D
LENS DESIGN TECHNOLOGY:	3-Zone Progressive™ Design, center-near aspheric optics
BASE CURVE:	8.5 mm
DIAMETER:	14.2 mm
CENTER THICKNESS:	0.07 mm @ -3.00D
SPHERICAL POWERS:	+6.00D to -10.00D in 0.25D steps (including plano)
ADD POWERS:	Low: +0.75D to +1.50D spectacle Add High: +1.75D to +2.50D spectacle Add
VISIBILITY TINT:	Light blue
MODALITY:	Monthly; Daily wear indication

See package insert for more information.

REFERENCES: 1. Data on file. Bausch & Lomb Incorporated, Rochester, NY; 2015. 2. Results of an online survey with patients that wore Bausch + Lomb ULTRA® for Presbyopia lenses for approximately 5 days (n=395). Survey questions were top 3-box scores (% Strongly Agree, Agree, Slightly Agree) on a 6-point agreement scale.

®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates.
©2017 Bausch & Lomb Incorporated. UFP.0059.USA.17

BAUSCH + LOMB
See better. Live better.

Bausch + Lomb ULTRA® for Presbyopia FITTING GUIDE

STEP 1: Update spectacle refraction and Add power

STEP 2: Select contact lens distance prescription based upon spherical equivalent from spectacle Rx and following Add guidance (adjusted for vertex distance if necessary)

ADD SELECTION:

SPECTACLE Add	BOTH EYES
+0.75D to +1.50D	Low Add
+1.75D to +2.50D	High Add

EVALUATE THE LENS FOR SUCCESS

- Allow trial lenses to equilibrate for at least 10 minutes before assessing fit and vision
- Evaluate distance and near vision binocularly in normal room illumination
- If vision at distance and near are satisfactory, dispense lenses and schedule follow-up exam within 1-2 weeks



REFINE IF NEEDED

Determine eye dominance at distance by placing a +1.50 loose handheld trial lens alternately over each eye binocularly through updated distance correction. The eye for which binocular vision is blurriest through the +1.50 is the dominant eye.

EASY, PREDICTABLE FITTING¹

92% of patients

agree that Bausch + Lomb ULTRA[®] for Presbyopia contact lenses allow them to shift focus naturally from near to far throughout the day²

IF PATIENT IS WEARING:

		NEAR VISION		DISTANCE VISION		
		DOMINANT EYE	NON-DOMINANT EYE	DOMINANT EYE		NON-DOMINANT EYE
TWO LOW ADDS	Initial Lens	Low Add	Low Add	Initial Lens	Low Add	Low Add
	Refinement 1	Low Add	High Add	Refinement 1	Bausch + Lomb ULTRA[®] sphere	Low Add
	Refinement 2: If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye (wearing High Add lens) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.	Refinement 2: If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (wearing Bausch + Lomb ULTRA [®] spherical lens) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.				
TWO HIGH ADDS		DOMINANT EYE	NON-DOMINANT EYE	DOMINANT EYE		NON-DOMINANT EYE
	Initial Lens	High Add	High Add	Initial Lens	High Add	High Add
	Refinement 1	High Add	Add +0.25D to the non-dominant eye	Refinement 1	Low Add	High Add
	Refinement 2: If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.	Refinement 2: If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (wearing Low Add lens) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.				